Original Copy

Control Copy 2

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(Standard Operating Procedure)

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Clinical Pathology Laboratory Workflow and Project Implementation

1 Purpose

This SOP describes the workflow of clinical pathology laboratories and standardizes the implementation of clinical pathology items for safety evaluation studies.

2 Scope

Applicable to all staff of clinical pathology laboratory.

3 Role

Role	Responsibility
Department Head	Ensure that the personnel participating in the trial have been appropriately trained and passed the assessment, and have a general grasp of the work of the department. Be responsible for coordinating the content communication with the study director in clinical pathology and the operation of the project to ensure the authenticity and validity of the test data of the project and the scientificity of the analysis report.
Clinical Pathology Technician	Check the sample status and identification in strict accordance with each SOP and check the number of samples according to the sampling record. Sample test, record arrangement, instrument maintenance and other operations to ensure the accuracy of test results and good instrument status.

4 Definitions

NA.

5 Procedures

5.1 Project Implementation Procedures

5.1.1 Pre-Project Preparation and Pre-Test Material Preparation

- (1) The personnel designated by the department head for the testing performed by clinical testing department shall participate in protocol training before the trial, and put forward suggestions on the content and process of clinical pathology-related testing.
- (2) After receiving the study protocol and study schedule by mail, develop the test plan according to the study protocol, and make preparations before the trial, such as method development and verification, and storage of reagents and consumables.
- (3) According to the actual situation of the study, decide whether it is necessary to explain or train the research background, drug characteristics, detection cycle, expected results, detection of special indicators of the trial and whether the samples need special preservation

after the completion of detection.

(4) For the study using Provantis system to collect the clinical test data (except for screening animals during quarantine period, veterinary diagnosis and treatment and additional testing required by SD and other tests defined in special protocol), an application needs to be submitted in Provantis. The technician logs in to Provantis system, enters "clinical pathology" module, selects "application for scheduled test" or "application for unscheduled test", selects the study with test application, date, animal, item, etc., to generate the test application. When submitting an application in Provantis, the protocol definition in Provantis must be performed in accordance with the items specified in the protocol, and do not over-choose or omit the items to ensure that the application information is correct.

(5) Prepare sample labels, collection tubes, and related record forms in advance based on the testing application information. When preparing the collection tube by pasting the label, pay attention to checking the sample information on the label (topic number, animal number, date, sample type) and the quality of the collection tube (sealing, whether there is damage, amount of anticoagulant, etc.) to ensure that the information is correct and the quality is qualified. If necessary, according to the collection volume requirements and anticoagulant ratio requirements of the plan, prepare anticoagulant tubes (such as sodium citrate anticoagulant tubes) and fill out the "Reagent Preparation Record Form (SPAT-011-F1)".

(6) Verification of applied testing indicators: After applying for testing indicators in the Provantis system, a second person needs to verify the testing indicators and animal information according to the experimental plan to ensure accuracy, and fill out the QC form (SMAN-044-F1).

5.2 Receipt and Handling of Specimens

5.2.1 Specimen Receipt

Upon completion of specimen collection, the Clinical Pathology Laboratory Technician is responsible for receiving the specimen. Refer to SMAN-033 "Study Materials Transport And Reception Record" for details. The sample receiving personnel shall check the sample test number and animal number, sample size and sample status (the sample collection volume shall meet the requirements specified in the protocol, there is no clot in anticoagulated whole blood, and the serum/plasma shall be light yellow/yellow transparent/translucent liquid after centrifugation). Following completion of reconciliation, appropriate documentation will be documented on the "Pre treatment record of clinical test samples (SCLI-001-F1)".

If unqualified samples are found (for example, the information on the sample tube is inconsistent with that on the application information; the collection container is incorrect, for example, the anticoagulant tube is used for serum samples; the samples are not submitted in time according to SOP requirements, resulting in failure to process and test the samples within the specified time; the sample size is insufficient by visual inspection or the sample suction volume is insufficient by the instrument during on-board testing, resulting in failure to detect the samples normally; the anticoagulant sample is observed by visual inspection to have clots or the sample is indicated by the instrument during on-board testing, resulting in abnormal

results or curves; the serum/plasma is found to have severe hemolysis, jaundice or lipemia after centrifugation, or the corresponding information is indicated by the instrument during on-board testing, resulting in abnormal results; etc.), the clinical technicians shall timely contact the study director (SD). If re-collection is required, the SD shall notify the test technician responsible for the study after approval. Recollected specimens were noted as recollected in Provantis or on the "Clinical Pathology Specimen Pretreatment Record (SCLI-001-F1)". Where some unsuitable samples still need to be tested after consultation with SD due to some special reasons such as difficulty in obtaining samples, the unqualified status of samples shall be indicated in the Results Remarks Information Column of Provantis System. Rejected sample information is documented on the "Pre treatment record of clinical test samples" (SCLI-001-F1).

5.2.2 Specimen Processing

Refer to SCLI-004 "Preparation and numbering of clinical pathological specimens" for the processing of different types of specimens

5.3 Specimen Testing

- (1) Check the test indicators: before sample test, the operator checks whether the indicators are correct according to the study protocol.
- (2) Instrument preparation: Perform startup maintenance of the instrument before the start of the test and record it on the relevant instrument use record form. Ensure that the instrument is prepared before sample analysis.
- (3) Calibration and Adjusting: For quantitative test items such as biochemical analyzer, coagulometer and hemocytometer, it is required to re-calibrate or calibrate after the reagent changes the batch number or after the instrument maintenance replaces the key components. The data generated by calibration shall be printed and signed for confirmation in real time, and archived annually with the instrument use records.
- (4) Control detection: Control shall be used in strict accordance with the requirements in the instructions for use. In order to ensure the quality of testing, quality control of corresponding instruments should be performed every day before the start of sample testing. If there are more than 20 samples detected by the corresponding instrument on that day, or the instrument has been running for more than 8 hours, an additional quality control of shutdown should be conducted after the end of the experimental testing. Two or more concentration levels should be selected for quality controls (except for special items or quality controls whose corresponding concentration cannot be provided by the manufacturer). Samples can be detected only when the quality control results are qualified before the start of sample detection, and the shutdown quality control result (if any) must be qualified before the test of the day can be ended. When the doubtful circumstances (such as instrument malfunctions, reagent problems, etc.) are noted, additional quality control test shall be performed. If the quality control before or after the sample test on that day fails, the cause of loss of control shall be analyzed and the test shall be repeated. After correcting the quality control loss during or after the experiment, it is necessary to evaluate and process the samples tested before the

loss of control, and report to the laboratory manager for review. Reasons for quality control failure and corrective actions for quality control need to be described in detail in the Provantis system, which also should be documented on the "Record of Loss of Control of Internal Quality Control (SCLI-001-F2)" if instruments unable to use the Provantis system. According to SCLI-002 "Clinical Pathology Quality Control Serum Reconstitution Preparation and Internal Quality Control Rules", quality control analysis, evaluation and out-of-control handling should be performed.

(5) Testing: Perform testing and analysis according to the SOP of each instrument.

When samples are detected on the machine, the sample barcode is automatically identified by the instrument and detected according to the identification order. If the instrument cannot identify the sample barcode and manual entry is required, the instrument analysis number is prepared by the clinical technician. For each test, the corresponding situation list of animal number and instrument analysis number should be described and recorded in the "Record of Pretreatment of Clinical Pathological Specimens" (SCLI-001-F1). Performing the test. If the Study Director has specific requirements for the order of specimen testing, this should be stated in the protocol.

5.4 Review of Results

The inspection related personnel shall review the results after the completion of detection. Carefully check whether the test results have such conditions as exceeding the test line, sample status alarm and abnormal instrument status alarm. Specimens with abnormal results were retested according to SCLI-003 Criteria for Repeated Determination of Clinical Pathology Specimens. Technicians shall timely report abnormal results to the study director or department head when discovering them, and timely report abnormal results with typical characteristics in the test to the study director.

5.5 Release of Results

- (1) For trials not involving Provantis system, a paper report will be printed after verifying the results on the instrument. After the operator and reviewer have signed and confirmed the paper report, it will be submitted as the original data to the study director.
- (2) For the collection of clinical pathological data using Provantis system, the detection results of each instrument will be automatically transmitted in real-time to the Provantis system. When real-time transmission is not possible due to special reasons, a reviewer needs to review the results in the instrument to ensure that they are correct before manually transmitting or recording them back to the Provantis system. After the test personnel comment on the results, data quality identification and data substitution identification, the reviewer reviews the results and electronically signs them to ensure that the results are correct and then publishes the results. After the results are released, the study director can consult the clinical test data in the Provantis system test module.

5.6 Sample Storage

Serum or plasma should be stored according to the protocol. If dispensing is required for

elective testing, it should be recorded in the “Clinical Pathological Specimen Sub-packing Record (SCLI-001-F3)”. The storage period should be at least until the study is completed to determine the specimen processing with the sponsor. Other types of specimens may not be preserved if specimen preservation is not specified in the study. Specimen storage and processing shall be documented on the “Sample Keeping Records of Clinical pathology (SCLI-001-F4)” and “Sample Handling Records of Clinical pathology (SCLI-001-F5)”, respectively.

5.7 Records Collation and Archiving of General Information

(1) After the completion of the test record, the tester shall sign and confirm it, and check it by the checker and sign it. On that day, the records are sorted out and temporarily stored by the testing technician of Clinical Pathology Department or delivered to the study director for storage.

(2) General data of clinical pathology department include original study records, reagent reconstitution records, specimen preservation records, specimen processing records, quality control out-of-control records, use and maintenance records of detection instruments, instructions of test reagents (the instructions of the same version are only archived once), delivery test reports of reagents (the same kind of reagents are only archived once), instructions and target values of controls and standards, SOP deviation records, inter-laboratory quality assessment (capacity verification) records and unsatisfactory results processing records, and other records generated in the study. General data of clinical pathology department were archived annually, and original records were archived as per the study.

6 Handling of Unexpected Events

A. Emergency handling of instrument malfunctions

Urgently contact the instrument engineer for remote guidance or on-site inspection and maintenance;

Sample testing: After interruption of testing, store the samples properly (according to SOP "SCLI-004-V3.0 Clinical Pathology Standards Sample storage requirements in the preparation and numbering of this document; If the instrument cannot be restored to normal for a long time, the sample should be promptly sent to a backup cooperative laboratory for testing;

Recovery testing: After repairing the instrument malfunction, it is necessary to calibrate and quality control the instrument to ensure that its performance is normal before resuming sample testing. After major fault repairs, engineers evaluate whether instrument performance validation is necessary, and only resume sample testing after validation is passed.

B. Emergency response to power and water outages

Assess the duration of power and water outages: Contact the operations department to confirm the duration of power outages, short-term power outages, and UPS power supply

Startup can provide 30-60 minutes of buffer power supply for the detection instrument; The laboratory needs to keep a certain amount of distilled water/pure water on hand for emergency preparation of reagents and cleaning;

Protecting samples and data: When the power outage exceeds the support time of the UPS

power supply, sample testing should be paused, and the running instruments should be urgently shut down according to SOP to maximize instrument protection and data preservation. All interrupted sample testing batches should be recorded;

Sample testing: After the testing is interrupted, store the samples properly (according to the sample storage requirements in SOP "SCLI-004-V3.0 Preparation and Numbering of Clinical Pathological Specimens"); If the instrument cannot be restored to normal for a long time, the sample should be promptly sent to a backup cooperative laboratory for testing;

Notification and Communication: Immediately notify the laboratory manager, SD, etc. to assess the impact on the project timeline;

Recovery testing: After the water and electricity are restored, start the equipment in sequence, execute the startup program, perform startup maintenance, and perform quality control testing to ensure that the instrument is in normal condition and the quality control is qualified before resuming sample testing.

7 Relevant Documents

Document Number	SCLI-004	Document Title	Preparation and numbering of clinical pathological specimens
	SCLI-002		Clinical pathology quality control serum reconstitution preparation and internal quality control rules
	SCLI-003		Criteria for Repeated Testing of Clinical Pathology Specimens
	SMAN-033		Transfer and Receipt of Test Materials
Form Number	SCLI-001-F1	Form Title	Retreatment record of clinical pathological specimens
	SCLI-001-F2		Indoor QC runaway record
	SCLI-001-F3		Clinical Pathology Sample Dispensing Record
	SCLI-001-F4		Clinical Pathology Specimen Retention Record
	SCLI-001-F5		Clinical Pathology Specimen Processing Record
	SPAT-011-F1		Reagent Preparation Record Form
	SMAN-044-F1		Quality Control Form
Template Number	NA	Template Title	NA

8 References

- (1) The United States Food and Drug Administration (US FDA) Good Laboratory Practice for Nonclinical Laboratory Studies (21 CFR Part 58).
- (2) OECD Good Laboratory Practice for Nonclinical Laboratory Studies.
- (3) National Medical Products Administration (NMPA) Good Laboratory Practice for Nonclinical Laboratory Studies.

9 SOP Revision History

Revision History

SOP current version No.	Effective date of the old version	Reviser	Description of revision
V2.0	2024-05-14	Xiaoqing Ren	1. 5.3 In specimen testing, item 4 is modified to read: if the test number exceeds 40 or the first quality control time exceeds 8 hours, a quality control should be added after the test, and the operation can be finished after the quality control passes; 2. SOP unified bilingualism.
V3.0	2024-10-01	Xiaoqing Ren	1. Modified the content of quality control testing in section (4) of 5.3. 2. Added and modified the content of (4) and (5) in section 5.1.1. 3. Modified 5.5 release of results. 4. Added form in 6 Relevant Documents.
V4.0	2025-04-17	Xiaoqing Ren	1. Added the content of "verifying sample labels and sampling tube information" in section 1.5.1.1 (5). 2. Added (6) Verification of applied testing indicators in section 5.1.1. 3. Added 6 Handling of Unexpected Events. 4. Added Quality Control Form (SMAN-044-F1) in 7 Relevant Documents.