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1 目的	3
2 适用范围	3
3 职责	3
4 定义	3
5 操作程序	3
5.1 操作方法	3
5.2 拓扑图	4
5.3 账户管理	4
5.4 数据保存和备份	5
5.5 数据恢复	5
5.6 注意事项	5
5.7 维护保养	5
5.8 性能验证	5
5.9 故障处理	5
6 相关文件	5
7 参考文献	6
8 SOP 修订历史	6
1 Purpose	7
2 Scope	7
3 Role	7
4 Definitions	7
5 Procedures	7
5.1 Operating Procedure	7
5.2 Topological map	8
5.3 Account management	8
5.4 Data retention and backup	9
5.5 Data recovery	9
5.6 Notes	9
5.7 Maintenance	9
5.8 Performance Qualification	9
5.9 Troubleshooting	10

6 Relevant Documents	10
7 References	10
8 SOP Revision History	10

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Use and Maintenance of Cellometer K2 Fluorescent Cell Counter

1 Purpose

The purpose of this SOP is to standardize the methods for using, maintaining, and troubleshooting of the Cellometer K2 Fluorescent Cell Counter.

2 Scope

This SOP is applicable to all personnel using the Cellometer K2 Fluorescent Cell Counter at this facility.

3 Role

Role	Responsibility
Operators	Must understand and comply with this SOP.
Equipment Manager	Responsible for operating and maintaining the equipment according to this SOP.
Administrator	Responsible for creating accounts and assigning user rights, creating projects, assigning project-related personnel.

4 Definitions

NA.

5 Procedures

5.1 Operating Procedure

- (1) Connect the USB data cable to the computer, plug in the power, turn on the switch at the back of the equipment. The front indicator light of the equipment will illuminate accompanied by a sound.
- (2) Double-click the Cellometer Image Cytometer software icon on the computer desktop to launch the program and open the software interface.
- (3) Peel off the protective films from both sides of the counting chamber, add about 20 μ L of the sample to be measured, and insert the slide needing counting into the sample port.
- (4) Select the Assay Type based on the actual assay type:
For example: For Car-T (or T cells) dual fluorescence AO/PI counting, click on "Assay" in the "SETUP" tab, choose "Car-T" from the drop-down menu.

Sample ID

new sample

(5) Enter the sample name.

Preview B1 ☒ B1 ☐ B2

Preview F1

(6) Click "Preview Brightfield Image".

(7) Rotate the knob at the bottom right of the equipment to adjust focus.

(8) Cells should appear centered and bright, with clear outlines.

(9) Click the counting button  below to automatically capture images and count.

(10) Display the analysis results.

Results:

Count	Concentration	Mean Diameter
=====	=====	=====
Total Cells: 2412	5.25×10^7 cells/mL	3.5 micron
Live: 1852	4.03×10^7 cells/mL	3.6 microns
Dead: 560	1.21×10^7 cells/mL	3.1 micron


Viability: 76.9%

(11) Counting the next sample: Remove the counting chamber. If another sample needs to be counted, repeat Steps 3-10.

(12) Shutdown: Click the  in the upper right corner of the main page to exit the program, turn off the switch at the back of the equipment, and complete the Equipment Use and Maintenance Record (SMAN-024-F6).

5.2 Topological map

The network topology of this system is as follows:

5.3 Account management

The equipment uses a local account to log in, and account types include administrator and operator accounts. Each user is assigned a unique user name and set a password by the administrator, and uses his or her own user account to log in to the instrument to operate and fill in the equipment's

usage records. The password for the user account must be changed once at least every three months and the last three passwords cannot be reused. If no operation occurs within 5 minutes, the system will automatically lock.

(1) Administrator: The administrator account is usually held by a member of facilities operation department and is responsible for creating accounts and assigning user rights, etc.

(2) Operator: The operator account only allows operations such as editing and saving results.

5.4 Data retention and backup

(1) The data save path is: D:\K2 Data.

(2) Perform full NBU backups of the equipment's electronic data at least every 10 days.

(3) The original records of this equipment are paper records.

5.5 Data recovery

When data corruption requires data recovery, you should contact the administrator to make a backup of the current database, then copy the database files to be recovered to the designated save location, recover the corrupted files, and make a record.

5.6 Notes

(1) Before using, ensure the USB data cable is properly connected, the equipment is powered on, and open the software again.

(2) Remove both protective films from the counting chamber to prevent a dirty background affecting counting results.

(3) Before counting, thoroughly mix the sample; it is recommended to pipette and mix approximately 10 times to ensure accurate counting.

(4) PI dye can be counted immediately without reaction time after mixing with the sample.

(5) When aspirating samples or dyes, ensure to change the pipette tips; samples and dyes should be thoroughly mixed in the cell chamber.

5.7 Maintenance

(1) Clean the external surface of the cell counter after each use.

(2) If not used for an extended period, power on the equipment at least every six months to check system status, software, and hardware functionality. Record and complete the Equipment Use and Maintenance Record (SMAN-024-F6) accordingly.

5.8 Performance Qualification

Perform the performance qualification (PQ) annually to validate equipment performance. Ensure validation data is stored and archived as required.

5.9 Troubleshooting

In case of damage to accessories or system failure, the system should be managed and repaired with reference to SOP SMAN-024 Equipment Management and SMAN-062 Computerized System Emergency Management, and fill in SMAN-024-F8 Equipment Repair Record and SMAN-062-F1 Computerized System Exception Reporting and Processing Record.

6 Relevant Documents

Document Number	SMAN-024	Document Title	Management of Equipment
Form Number	SMAN-024-F6	Form Title	Equipment Use and Maintenance Record
	SMAN-062-F1		Computerized System Exception Reporting and Processing Record
	SMAN-024-F8		Equipment Repair Record
Template Number	NA	Template Title	NA

7 References

- (1) The United States Food and Drug Administration (US FDA) Good Laboratory Practice for Nonclinical Laboratory Studies (21 CFR Part 58).
- (2) OECD Good Laboratory Practice for Nonclinical Laboratory Studies.
- (3) National Medical Products Administration (NMPA) Good Laboratory Practice for Nonclinical Laboratory Studies.

8 SOP Revision History

Revision History

SOP current version No.	Effective date of the old version	Reviser	Description of revision
V2.0	2024-07-09	Yuwei Li	1. Bilingual SOP 2. “5.2 Notes” Add (6) Electronic data from cell counts need to be archived with the test.
V3.0	2024-10-17	Yuwei Li	1. Add 3 administrator of role, Add 5.2 Account management, 5.3 Project creation and audit trail, 5.4 Data retention and backup, 5.5 Data recovery
V4.0	2024-11-21	Yuwei Li	1. Changed the name of the equipment to Cellometer K2 Fluorescent Cell Counter. 2. 5.4 Data retention and backup time

			modified to at least once a month. 3. 5.4 Add “the original records of this equipment are paper records”.
V5.0	2025-02-27	Yuwei Li	1. Delete “5.3 Project audit trail”
V6.0	2025-04-24	Yuwei Li	1. Add “5.2 Topological map”. 2. Modify the backup of electronic data of the equipment in “5.4 Data retention and backup”. 3. modify the handling method of “5.9 Troubleshooting”.
V7.0	2025-06-20	Yuwei Li	1. Delete “5.3 Account management “(1) Administrator: creating projects, assigning project-related personnel and “5.6 Notes” (6) Electronic data from cell counts need to be archived with the test.